
NanoProtect®
SilverMask

PROTECTIVE MASK WITH FILTER

WITH SILVER ION

4-layer anti-bacterial and anti-fungal protective mask with a pocket for replaceable hydrophobic nanosilver filters

Our LionMED NanoProtect SilverMask protective mask with a pocket for hydrophobic filters with active nanosilver particles is a carefully made of polyester yarn with the addition of elastane, a perfectly fitting to the anatomical face shape, a protective mask tightly covering the upper roads Respiratory minimizing the risk of infection thanks to the antimicrobial properties of the LionMED filter exchange nanosilver. Pocket allows use of the attachedexchangeable two-layer filter made of non-woven fabric protected with bactericide and hydrophobic agent by using, among others silver ions that limit the multiplication of the microbiome captured from the air during filtration.

In the technological process, the filters used undergo sterilization process at 200 degrees Celsius, and the reusable mask itself is adapted for steam sterilization or washing at 95 degrees Celsius, thus allowing the mask to be kept at the highest purity.

The most important features of the mask:

- ✚ Reusable
- ✚ Can be sterilized at 65 C - 95 C
- ✚ 4-layer protection (when using a replaceable filter)
- ✚ Ultralight and comfortable even after long use
- ✚ Effectively inhibits the inhalation of bioaerosol
- ✚ Composition: 95% Polyester, 5% Elastane
- ✚ Size: universal, highly elastic
- ✚ Possibility of using replaceable filters
- ✚ Low cost of use 1 day - 1filter
- ✚ OEKO-TEX certificate of raw material

The most important features of the filters:

- ✚ For single use, max. 8 hours in a normal environment
- ✚ Neutralizes hostile microorganisms such as bacteria and fungi
- ✚ It blocks the uncontrolled multiplication of viral cells
- ✚ Ionized nanosilver - effectively reduces the risk of developing pathogenic flora
- ✚ 2 layers
- ✚ Size: 100x140mm filter area
- ✚ OEKO-TEX certified nonwoven filter

**PROTECTIVE
MASK WITH FILTER**

WITH SILVER ION


EU
manufactured

NanoProtect®
SilverMask

PROTECTIVE MASK WITH FILTER

WITH SILVER ION

NanoProtect SilverMask filter cartridge

Is a nonwoven fabric prepared with a view to maintaining purity, hygiene and resistance to the bacteria that surround us. The nonwoven fabric was protected with a specialized bactericidal and hydrophobic agent. It contains silver ions that control bacteria with a triple mechanism.


The antibacterial finish of the nonwoven with a silver-based agent prevents unrestrained growth of bacteria on its surface. Silver nanoparticles exhibit bactericidal and fungicidal activity. The natural balance of bacterial flora on the skin of the body is not affected. The finish used is a non-migratory antibacterial agent. This means that it works on the product without creating inhibitory areas (according to the Agar-Diffusionsplattentest test (SN 195 920) To increase the volume of safety against droplet-borne bacteria, the non-woven fabric is equipped with an additional hydrophobic coating that prevents liquid from penetrating inside the non-woven structure. In the presence of moisture, they are released antibacterial silver ions, which inhibit further multiplication of bacterial cells, protecting our health.

How is it working?

Silver nanoparticles are able to kill microorganisms, which is why they are currently used in the production of modern disinfectants. Nanosilver is an effective destructive agent against a wide spectrum of Gram-negative and Gram-positive bacteria, not excluding antibiotic-resistant strains. Nanosilver has a strong positive charge, which means it can literally search for pathogens. This is a property that the antibiotics used so far did not have. Silver nanoparticles do not have to rely on the case that they just come into contact with a cell of a pathogenic organism. It immediately destroys an enemy cell through oxidation.

Description of NanoProtect SilverMask use:

- 1 Wash your hands thoroughly before use
- 2 Take the mask and the attached filter from the packaging
- 3 Slide the unfolded filter into the pocket of the mask so that it does not protrude beyond its edges
- 4 Stretch the mask handles and put them behind the ears pressing the mask against the face so that it covers the face tightly over its entire surface
- 5 If necessary, adjust the filter in the mask pocket so that it covers the mouth and nose

Replacement of filters and safety rules

- + Each time the filter is changed, it must be accompanied by sterilization of the reusable protective mask. Never re-use the mask without sterility and safety.
- + Always dispose the used filter safely for yourself and the environment using disposable gloves.
- + Remember that the filter is disposable and attempting to reuse it may endanger your health and life.

Mask maintenance

- + Wash or sterilize at 65 - 95 ° C
- + Do not bleach and chlorinate
- + Do not dry mechanically
- + Do not iron

95% PES 5%LYCRA

PROTECTIVE
MASK WITH FILTER

WITH SILVER ION


EU
manufactured

NanoProtect[®]
SilverMask

PROTECTIVE MASK WITH FILTER

WITH SILVER ION

We provide full product photography and sales support

+
PROTECTIVE
MASK WITH FILTER

WITH SILVER ION


EU
manufactured

NanoProtect[®]
SilverMask

PROTECTIVE MASK WITH FILTER

WITH SILVER ION

**4-layer anti-bacterial and anti-fungal protective mask
with a pocket for replaceable hydrophobic nanosilver filters**

in 1 package
Mask + filter

Packing

Name	Quantity in unit package	Packaging size.(mm)	collective pack.	collective pack. size (cm)
4-layer anti-bacterial and anti-fungal protective mask with a pocket for replaceable hydrophobic nanosilver filters	1 pcs.	100x80x10	50 pcs.	30x22,5x10

Palletizing

Type of pallet	pcb quant. on the pallet	nb. of collective pack. on a pallet
120x80 EUR	9600 pcs.	192 cartons

**PROTECTIVE
MASK WITH FILTER**

WITH SILVER ION

Protective face mask
MO.SILVER.OP1

1 pcs.


EU
manufactured

NanoProtect[®]
SilverMask

PROTECTIVE MASK WITH FILTER

WITH SILVER ION

Replaceable LionMED NanoProtect SilverMask Filter available in
10 pcb separate packs

in 1 package
10 pcs. filter

Packing

Name	Quantity in unit package	Packaging size.(mm)	collective pack.	collective pack. size (cm)
LionMED filter for NanoProtect SilverMask 10-pack	10pcs.	140x100x10	50 sets.	30x29x12

Palletizing

Type of pallet	pcb quant. on the pallet	nb. of collective pack. on a pallet
120x80 EUR	6200 pcs.	124 cartons

PROTECTIVE
MASK WITH FILTER

WITH SILVER ION

LionMED filter for
NanoProtect SilverMask

10pcs.


EU
manufactured

NanoProtect®
SilverMask

PROTECTIVE MASK WITH FILTER

WITH SILVER ION

OEKO-TEX certificate for yarn used for the production of the NanoProtect mask

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
PROTECTIVE MASK WITH FILTER

WITH SILVER ION

OEKO-TEX Passport for pigment inks used in the production of NanoProtect filters

Nissenken
Quality Evaluation Center
OEKO-TEX® Laboratory
2-15-11 Kurumae
Taisho-ku
113-0001 Tokyo, Japan

OEKO-TEX®
CONFIDENCE IN TEXTILES

CERTIFICATE

The Company

SEIKO EPSON CORPORATION
80 Harashinden, Hirooka, Shiojiri-city
JP-Nagano, JAPAN

is granted authorisation according to ECO PASSPORT by OEKO-TEX® to use
the OEKO-TEX® mark

OEKO-TEX®
CONFIDENCE IN TEXTILES
ECO PASSPORT

NEP 1601 Nissenken

Textile chemicals. Tested and verified.
www.oeko-tex.com/ecopass

for the following chemical products

Product(s): See attached enclosure

Category: Disperse dyes

Supporting documents

- Declaration of conformity in accordance with EN ISO 17050-1 included in ECO PASSPORT by OEKO-TEX® Terms of Use.
- Analytical test report number: NEP 1601-2
- RSL Screening Report
- Detailed information about the components and safety data sheets of the chemical products mentioned above.

The above captioned product(s) can be used for the production of human-ecological optimized textiles. The combined results of the reports mentioned above reveal that there is no harmful effect on the human and environmental health of the textiles treated/finished with the above mentioned products. This evaluation used the test methods and requirements of the ECO PASSPORT by OEKO-TEX® that were in force at the time of the evaluation date. ZDHC MRL Conformance Level 1 is achieved for certified product(s) without restriction(s).

The certificate NEP 1601 is valid until 31.10.2019

Japan, 30.10.2018

Kazuhiro Funakawa
General Manager, OEKO-TEX® Laboratory
Nissenken

OEKO-TEX® Association | Fensterstrasse 23 | P.O. Box 2006 | CH-8027 Zurich

113

**PROTECTIVE
MASK WITH FILTER**

WITH SILVER ION


EU
manufactured

NanoProtect[®]
SilverMask

PROTECTIVE MASK WITH FILTER

WITH SILVER ION

OEKO-TEX Passport for pigment inks used in the production of NanoProtect filters

Nisshinkan
Quality Evaluation Center
OEKO-TEX® Laboratory
2-18-11 Kuramae
Taito-ku
111-0051 Tokyo, Japan
Company

OEKO-TEX®
CONFIDENCE IN TEXTILES

SEIKO EPSON CORPORATION
80 Harashinden, Hirooka, Shiojiri-city
JP-Nagano, JAPAN

ECO PASSPORT by OEKO-TEX® Lab Report No.

NEP 1601-2

ECO PASSPORT by OEKO-TEX® Certified Products

No.	Product name	Trade name	*Restriction(s)	ZDHC LEVEL
1	Ultrachrome DS Ink (Magenta), MAGENTA GENESTA® SBL 3L, MAGENTA GENESTA® SBL 10L	Ultrachrome DS Ink (Magenta), MAGENTA GENESTA® SBL 3L, MAGENTA GENESTA® SBL 10L	none	Level 1
2	Ultrachrome DS Ink (Black)	Ultrachrome DS Ink (Black)	none	Level 1
3	Ultrachrome DS Ink (Cyan), CYAN GENESTA® SBL 3L, CYAN GENESTA® SBL 10L	Ultrachrome DS Ink (Cyan), CYAN GENESTA® SBL 3L, CYAN GENESTA® SBL 10L	none	Level 1
4	Ultrachrome DS Ink (High Density Black), BLACK GENESTA® SBL 3L, BLACK GENESTA® SBL 10L	Ultrachrome DS Ink (High Density Black), BLACK GENESTA® SBL 3L, BLACK GENESTA® SBL 10L	none	Level 1
5	Ultrachrome DS Ink (Yellow), YELLOW GENESTA® SBL 3L, YELLOW GENESTA® SBL 10L	Ultrachrome DS Ink (Yellow), YELLOW GENESTA® SBL 3L, YELLOW GENESTA® SBL 10L	none	Level 1

Issue Date 30.10.2018

*Restriction(s): The parameter(s) mentioned under Restriction(s) have to be checked on the treated textile for compliance with the regulations of STANDARD 100 by OEKO-TEX®.

OEKO-TEX® Association | Genfstrasse 22 | P.O. Box 2006 | CH-8027 Zurich
Phone: +41 44 501 26 00 | info@oeko-tex.com | www.oeko-tex.com

3/3

**PROTECTIVE
MASK WITH FILTER**

WITH SILVER ION


EU
manufactured

NanoProtect®
SilverMask

PROTECTIVE
MASK WITH FILTER

WITH SILVER ION

**Certification and technical specification
of ®RUCO-BAC AGP nanosilver ionizing
filter LionMED**

+
PROTECTIVE
MASK WITH FILTER

WITH SILVER ION


EU
manufactured

NanoProtect[®]
SilverMask

PROTECTIVE MASK WITH FILTER

WITH SILVER ION

Certification of the ®RUCO-BAC AGP silver ionizing LionMED antibacterial filter

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